Bedienungsanleitung
Instruction manual
Manual de instrucciones
Mode d'emploi

Eltego


DE

EN

ES

FR

Enjoy mobility.

Inhalt

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[illegible]

Content

1.	Introduction	46
1.1.	Signs and Symbols	47
2.	Important Safety Instructions	47
2.1.	General Safety Instructions	47
2.2.	Safety in Vehicles	49
2.3.	Advice Regarding EMC Interference	49
2.4.	Brakes	50
2.5.	Use on Public Roads	51
2.6.	Purpose	51
2.7.	Indications for Use	52
2.8.	Contraindication	52
2.9.	Declaration of Conformity	52
3.	Description of Product and Content	53
3.1.	Inspection on Delivery	53
3.2.	Content	54
3.3.	Summary	54
3.4.	Type Plate and Serial Number	55
3.5.	Basic Settings	56
3.5.1.	Leg Supports	56
3.5.2.	Backrest	57
3.5.3.	Backrest Sag	57
3.5.4.	Suspension	58
3.5.5.	Armrests	58
3.5.6.	Side Parts	59
3.5.7.	Control Unit	59
3.5.8.	Lap belt	60
4.	Using the Wheelchair	61
4.1.	Description of Control Functions	61
4.2.	Getting In and Out	63
4.3.	Driving the Wheelchair	64
4.4.	Push Mode Operation	66
4.5.	Leg Supports	67
4.6.	Anti-tipping Wheels	67
4.7.	Lighting System	68
4.8.	Rear Reflectors	68

5.	Acessories	69
5.1.	Head Rest	69
5.2.	Crutch Holder	69
6.	Technical Data	70
7.	Maintenance Tips for Users	72
7.1.	Cleaning and Disinfection	72
7.2.	Check Before Use	73
7.2.1.	Brake Inspection	73
7.2.2.	Checking the Lights and Horn	74
7.2.3.	Inspection of Tires	74
7.3.	Batteries	74
7.3.1.	Charging Procedure	75
7.3.2.	Exchanging Batteries	76
8.	Maintenance Instructions	77
9.	Maintenance Tips for Retailers	78
10.	Parking	78
11.	Disposal	78
12.	Use by Others	79
13.	Transport	79
14.	Storage and Shipping	79
15.	Service life	80
16.	Warranty	80

1. Introduction

Dear user,
You have decided on a high-quality Bischoff & Bischoff electric wheelchair, and we would like to thank you for your trust.

The electric wheelchair was designed to offer you many advantages and meet your needs. The wheelchair can be fitted to your individual body size thanks to a multitude of adjustable features and optimal adaptability such as the length of the lower leg, the recline of the backrest, the sag of the backrest as well as the swing away function of the control unit and the choice of mounting the control unit either to the left or the right side.

The electric wheelchair is designed for indoor and outdoor use. Before using your new electric wheelchair for the first time, please read and follow the instructions in this manual carefully. The instruction manual is an integral and necessary part of the wheelchair. This is why you should keep it close by and include it when passing on the wheelchair to someone else.

Any repairs and adjustments require special technical training and may therefore only be carried out by authorized Bischoff & Bischoff mobility dealers.

Note!

In spite of careful research and editing, incorrect information may have crept into our operating manuals. Rights to technical changes reserved. Illustrations may include accessories. The current version of each instruction manual is ready for download at the Bischoff & Bischoff website:
www.bischoff-bischoff.com

1.1. Signs and Symbols

Warning! These are especially safety-relevant instructions. Follow the directions in the manual!

Important! This is especially useful information on a given subject.

2. Important Safety Instructions

2.1. General Safety Instructions

Before setting off, check that all attachments are properly fixed in place.

- To avoid falls and dangerous situations, you should first practice driving your new electric wheelchair on level, unobstructed ground. An attendant is recommended in this case.
- Familiarise yourself thoroughly with the braking and acceleration behaviour when driving in a straight line and into and out of curves.
- Please note that the shifting of weight through body movements or loading can increase the risk of tipping over to the side, to the front or back.
- The electric wheelchair must only be used to transport one person. Please respect the intended use (see Chapter 2.6).
- When setting up or adjusting the product, be aware of the risk of body parts being pinched.
- Never use the electric wheelchair under the influence of alcohol or other agents which affect your awareness or your physical and mental capacities.

- It is essential that you use ramps to get over obstacles. If it should be necessary to lift the Ejoy over an obstacle, then only do this without an occupant. Only use fixed frame parts to lift the chair.
- Avoid driving into obstacles without braking (furniture, steps or similar items).
- The control unit must always be switched off whenever you are getting in or out of the chair.
- When on roads, observe the road traffic regulations.
- The lighting and reflectors on the wheelchair must always be visible to other road users.
- Where possible, wear bright conspicuous clothes. You will then be seen better by other road users.
- Please note that upholstery parts which are subject to sunlight will heat up and could cause skin injury when touched. You should therefore cover these parts or protect the wheelchair from direct sunlight.
- The seating upholstery and pads' materials meet the flammability resistance requirements in accordance with EN 1021-2. Therefore do not expose the upholstery or pads or any other parts to ignition sources - such as, for example, cigarettes.
- Always put the seat and backrest into an upright position when you are driving the wheelchair.
- Use on slopes is only permitted when the seat and backrest are in an upright position.
- At low ambient temperatures, the metal parts of the wheelchair can get very cold, so please avoid skin contact if this is the case.
- If you should notice your electric wheelchair behaving unusually or if it is damaged, stop using it and contact your dealer.
- Reporting incidents:

If you, as the operator, user or their relatives, detect possible defects or functional limitations, please contact an authorised dealer immediately. Operators, users or their relatives should inform the dealer who supplied the dealer who provided you with the product of any suspected serious incidents which directly or indirectly have had

or could have one of the following and by which they are affected. You can also report this to the competent higher federal authority report

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of the health of the state of health of a patient, user or other persons,
- a serious risk to public health.

Please note that operating the wheelchair in push mode by uncoupling the drive mechanism will deactivate the electromagnetic brakes. The wheelchair should therefore only be used on level ground in this mode. Non-observance may lead to life-threatening situations.

2.2. Safety in Vehicles

The Ejoy is not designed to be used as a vehicle seat. Under no circumstances may it be used as such. It does not fulfil the requirements in compliance with ISO 7176-19. Violation can cause the most serious injuries resulting in death in the case of an accident (Fig. A).

2.3. Advice Regarding EMC Interference

Despite meeting all applicable EMC guidelines and standards it is possible that the wheelchair will experience interference from other electrical devices (e.g. cell phones) or interfere with these. Should you notice such interference, either turn off the other device or your electric wheelchair and move it outside the interference range.

2.4. Brakes

Operating brake:

When in motion the motors act as an operating brake. To brake the wheelchair slowly, move the steering and driving lever (joystick) slowly back to the neutral position (initial position).

For more forceful braking, release the steering and driving lever.

The braking distance is significantly longer on slopes than on flat ground. Bear this in mind when going down ramps and brake in time.

Make sure that you brake your electric wheelchair in good time. Particularly when approaching people or when driving downhill.

Parking brake/parking:

Switch the wheelchair off using the on/off switch (see Chapter 4.1).

If the drives are engaged, you cannot push or move the electric wheelchair (see Chapter 4.4).

Pushing mode:

In order to be able to push the wheelchair, pull the release levers on both sides of the wheelchair upwards (see Chapter 4.4), if the wheelchair is switched on, then the battery status display will flash. For pushing mode you can switch the wheelchair off.

Please note that push mode by disconnecting the electromagnetic brakes must only be used on flat ground. The electromechanical brake is then disabled. Non-observance means danger to life.

Emergency brake:

Press the on/off switch to initiate an emergency stop when driving.

This will lead to intense deceleration. This can lead to injuries and is therefore only to be used in emergencies.

An emergency stop must only be carried out in an emergency as it can lead to injuries.

2.5. Use on Public Roads

The wheelchair is designed for indoor and outdoor use.

Please keep in mind that you are using public roads and therefore have to observe all road traffic regulations.

Do not endanger other road users through inconsiderate driving, especially on walkways.

2.6. Purpose

The Eltego electric wheelchair is designed for indoor and outdoor use by adult persons with impaired mobility.

It is therefore sufficiently compact and agile for indoor use and also capable of handling medium distances and overcoming small obstacles outdoors. It is classified under category B according to the EN 12184 standard.

The maximum load capacity is 120 kg.

2.7. Indications for Use

Provision of this electric wheelchair is suitable for adults with significant through to fully-developed impairment of mobility/ walking with structural and/or functional damage of the lower extremities, e.g. through

- paralysis
- loss of limb
- defect/deformity of limb
- joint contracture
- joint damage (not in both arms)
- other diseases

The provision of an electric wheelchair is indicated if the use of a hand-driven wheelchair is no longer possible but the proper operation of an electric drive is possible.

2.8. Contraindication

The provision of an electric wheelchair is unsuitable for adult persons

- with severe balance disorders
- with impaired or insufficient vision
- with strongly impaired cognitive abilities

2.9. Declaration of Conformity

We Bischoff & Bischoff GmbH declare on our sole responsibility, that the wheelchair Ejoy complies with all requirements of MDR 2017/745 which are applicable.

3. Description of Product and Content

3.1. Inspection on Delivery

All Bischoff & Bischoff products undergo a proper final inspection in our company and are labeled with the CE mark.

The wheelchair is delivered in a special carton. Keep the carton after unpacking, if possible. It is useful for later storage of the product or for sending it back. The Eltego wheelchair comes disassembled in its component parts.

Assembly, adjustment to your body size and first instructions are taken care of by your qualified mobility dealer.

If the wheelchair is shipped by rail or road, the goods should be checked immediately for transport damages (carton) in the presence of the bearer. Should there be any damage, please contact our customer service immediately.

Check the contents for completeness and make sure nothing is damaged. Should you find any irregularities or damage, please contact our customer service.

3.2. Content

After receipt of the goods, please check the contents for completeness right away. The contents consist of:


- outer package
- electric wheelchair (ready for use incl. 2 batteries)
- charging unit
- instruction manual

3.3. Summary

3.4. Type Plate and Serial Number

The type plate and serial number are located under the seat. The values given here are for illustration purposes only. Please refer to the type plate on your wheelchair for the exact specifications.

- (A) Manufacturer logo
- (B) Manufacturer information
- (C) Medical device
- (D) CE mark
- (E) Vehicle fixing
- (F) Warning! - read the User manual
- (G) Important! - read the User manual
- (H) Model designation
- (I) Model number
- (J) max. user weight
- (K) Date of manufacture
- (L) Classification in accordance with DIN EN 12184
- (M) maximum speed
- (N) UDI number
- (O) Serial number

3.5. Basic Settings

It is necessary to adjust the wheelchair to your body size before first use. This is the only way to ensure a correct posture allowing long drives without tiring. Below you find a description of the basic settings of the wheelchair. You should carry out these adjustments with the help of your mobility dealer. You need a set of 4 mm, 5 mm and 6 mm allen keys, a screw driver and ring spanners, 10 and 13.

3.5.1. Leg Supports

Length of lower leg supports

- The lower leg supports can be adjusted from 470 to 550 mm.
- To do so, loosen the two allen screws (SW 4).
- Adjust to the desired height (Fig. 3).
- Tighten the screws again. Use Loctite 243 to lock screws in place.

Make sure to always use two screws for secure fastening!

Angle of footrest

- The angle of the footrests is adjustable in 25° steps.
- Loosen the screwing at the side of the footrest.
- Adjust to the desired angle (Fig. 4).
- Tighten the screw again. Use Loctite 243 to lock screw in place.

Make sure to set the angle only within the marked range!

3.5.2. Backrest

The angle of the backrest of your electric wheelchair has a manual recline function and can be easily adjusted to the wanted position by pulling the strap at the lower end of the backrest. You also have the option of folding the backrest completely down to the front.

- You can adjust the angle of the backrest from 90° to 120° in four steps.
- Pull the cable behind the back while holding the backrest with the other hand (Fig. 5).
- After letting go of the cable the bolts will lock into place.


Make sure that both bolts on the left and right are completely locked in place.

3.5.3. Backrest Sag

A well-adjusted backrest makes it easier to relax and sit for a long time and effectively reduces the risk of developing pressure sores. Make sure there is enough sag in the pelvic area.

- Remove the padding from the tension straps of the backrest.
- Loosen the velcro fasteners on the straps.
- The straps are adjusted one after the other from the bottom up.
- Slide the padding over the tension straps again (Fig. 6).

3.5.4. Suspension

You can adjust the suspension to your weight for a comfortable drive.

- You will need a special wrench to adjust the spring deflection. (Please contact the service department)
- Pushing the spring together results in a reduced spring deflection and a "harder" setting.
- Expanding the spring increases the spring deflection and guarantees a "softer" setting (Fig. 7).

Make sure that both suspensions on the right and left are adjusted the same.

3.5.5. Armrests

The armrests can be completely taken off to be able to get in and out of the chair from the side. In addition, the height of the armrests can be adjusted. They can also be adjusted sideways by adjusting the side parts (see side parts) (Fig. 8)

Removing the armrests

- To remove the armrest, loosen the wing bolt and pull out the armrest to the top.
- When installing the armrest, make sure it is inserted all the way and fastened again.

Adjusting the height

- First take the armrest out of the mounting.
- Remove the arresting bolt and screw it back into the armrest at the desired position.
- The lowest hole represents the maximum height of the armrest.

3.5.6. Side Parts

You can slide the side parts sideways.

- Loosen the fastening (SW 13) below the seat and slide the side part into the desired position.
- Tighten the fastening again (Fig. 9).

Make sure the fastening is firmly tightened.


3.5.7. Control Unit


The control unit can be mounted either on the right or left-hand side underneath the armrest each (Fig. 10).

- The control unit must be adjusted by a qualified mobility dealer to warrant safe functionality of the wheelchair.

Adjusting the height

- The distance to the armrest can be adjusted.
- Loosen the screw connection on the inner side of the side part (Fig. 11).
- Move the control unit into the desired position, then tighten the screw again.

Distance from the armrest

- You can adjust the position of the control unit to your forearm length.
- Loosen the hand screw on the outer side of the side part. (Fig. 12).
- Move the control unit into the desired position, then tighten the screw again.

Don't forget that the cables to the control unit have to be repositioned. This should be done by a qualified mobility dealer.

3.5.8. Lap belt

The seatbelt is attached to the side of the backrest (Fig. 13).

- Adjust the length of the seatbelt by positioning the buckle accordingly.
- Fasten the lap belt with the latch.
- Open the lap belt by pushing the red button at the latch.



While using the wheelchair, the lap belt must always be properly fastened and closed.

4. Using the Wheelchair


After adjusting the wheelchair to your personal needs, familiarize yourself with the control functions.

4.1. Description of Control Functions


Battery Status Display (Fig. 16)

The battery status display shows to what extent the batteries are charged. If the indicator is within the green range, you have enough battery power available. If the indicator is within the red range, only short drives are possible and the batteries should be recharged. (See charging procedure).

ON/OFF Button (Fig. 17)

Pushing this button will put your wheelchair into operational mode.

Speed Range Display (Fig. 18)

The speed range display shows you the pre-selected speed range. You can set the speed between level 1 and 5, with the corresponding number of LEDs lighting up. The speed range indirectly tells you the maximum speed possible at full deflection of the joystick, i.e. with level 5 you can reach the maximum speed possible, while level 1 will keep you within the lowest speed range.

Setting the Speed Range (Fig. 18)

These buttons allow you to increase (right) or decrease (left) the maximum speed possible. Within each range you control the speed with the joystick.

Horn (Fig. 19)

The horn allows you to attract attention in difficult situations.

Immobilizer System

To protect your wheelchair against unauthorized use or theft, park it in a locked room or activate the immobilizer system.

- Push the ON/OFF button until you hear a beep, then let go.
- Move the joystick forwards, then backwards and then back into the neutral position. You will hear a beep each time. Now the control is blocked.
- Activate the control again by pushing the ON/OFF button.
- Move the joystick forwards, then backwards and then back into the neutral position. You will hear a beep each time. Now the control is activated again.

4.2. Getting In and Out

Getting in and out, sideways transfer

- Move the wheelchair as close as possible parallel to the other seat.
- Should you want to transfer from or to another wheelchair or rehabilitation stroller, secure this one first by activating the handbrake.
- Secure the Eltego against rolling away unintentionally by pushing the ON/OFF button and thus turning off the control (all displays cleared).
- Take off the armrest/side part.
- Fold the footrests to the side or take off the leg supports.
- Now slide sideways onto the other seat.
- Please make sure you are sitting on the seat as far back as possible.
- Put the armrest and leg supports back into their original position.

Warning! Risk to Life!
The drives may only be disengaged on level ground.

Getting in and out from the front

- First fold up the footrests or take off the leg supports.
- Move the Eltego as close as possible to the other seat so that they face each other.
- Should you want to transfer from or to another wheelchair or rehabilitation stroller, secure this one first by activating the handbrake.
- Secure the Eltego against rolling away unintentionally by pushing the ON/OFF button and thus turning off the control (all displays cleared).
- Now slide to the other seat by twisting your body.
- Finally, put the leg supports back into their original position.

4.3. Driving the Wheelchair

Before embarking on any journey, check that the magnetic brake is properly locked.

Startup

Sit down in your Bischoff & Bischoff wheelchair and rest your legs on the footrests.

Make sure you are sitting on the seat as far back as possible. Also make sure that no loose items of clothing such as scarves, shirts and the like are hanging from the side of the wheelchair, since they could be caught in the wheels!

Make sure you are sitting comfortably and securely. By pushing the ON/OFF button the Eltego is put into operational mode.

You can now use all electrical drive control functions.

Drive Control Functions

Your qualified mobility dealer can set the values for speed, acceleration and delay within the scope of statutory specifications.


Select a low speed range at first and then increase it during your drive, if necessary.

You can start driving by pushing the joy-stick forward. With the joystick you can also increase the speed and control the direction of travel. Putting the joystick in the neutral position will engage the brakes and secure the wheelchair against rolling away.

Driving straight ahead: To drive straight ahead, carefully push the joystick forward. The further you push the joystick forward, the faster the wheelchair will go. If you let go of the joystick you will automatically slow down.

Driving backwards: To drive backwards, carefully pull back on the joystick. The further you pull back on it, the faster your wheelchair will go. If you let go of the joystick, you will automatically slow down.

Turning: To change direction, push the joystick to the left or right. This is also possible when the joystick is in a forward or backward position, that is when you are going forwards or backwards.

4.4. Push Mode Operation

If you want to push the wheelchair manually in certain situations, the brakes have to be disengaged. In the free-wheeling mode the wheelchair is easy to push. The drive of the Eltego is equipped with mini-switches that deactivate automatically when free wheel mode is activated. When the on/off button is pressed, the battery indicator blinks nine times; the wheelchair can no longer be operated.

Warning! The electromechanical brakes then work no more

- Move the release lever up into the "Free Wheel" position (right and left).
- A small jolt will let you know that the drive is now disengaged.
- To re-engage the drive, push the lever back into the "Drive" position. It will lock automatically.

Warning! Always have both drives in the same position. (Risk of accident)

4.5. Leg Supports

To make it as easy as possible to get in and out of your wheelchair, you can fold the footrests up or pull the leg supports to the side and/or take them out.

- Open the locking mechanism to pull the leg support to the side. (Fig. 21)
- In the "outside position" the leg support can be taken off completely. (Fig. 22)
- To reinstall the leg support, repeat the steps in reverse order.
- Push the leg support back to the front until it catches with an audible click.

4.6. Anti-tipping Wheels

The anti-tilt rollers prevent the wheelchair from tipping over backwards.

- The anti-tipping wheels are especially recommended for novice or leg-amputated wheelchair users.
- The anti-tipping wheels are attached to the left and right side of the frame as a pair.

4.7. Lighting System

The electric wheelchair is equipped with a lighting system in accordance with the German traffic regulations (StVZO) for safe participation in street traffic (fig. 24). The lighting unit consists of 2 front reflectors, 2 side reflectors and 2 rear lights, which can also serve as reflectors simultaneously.

The active components are controlled with a button next to the joystick (fig. 25). Press the button briefly to switch the lighting system on/off. Hold the button down to adjust the brightness. The lighting system is independent of the controller and can also be operated when the wheelchair is switched off. When exiting the wheelchair, make sure that the light is switched off.

4.8. Rear Reflectors

The electric wheelchair is equipped with two rear reflectors which at the same time also serve as reflectors (Fig. 26).

Switch the rear on as soon as it starts to get dark, or when there are poor light conditions, so as to get a better perception of your surroundings.

They are attached by a clip to the frame at the back and can simply be pulled away to the side. By pressing the red button, you can switch on the rear reflex reflector, choose between different flashing modes and also switch it off again.

5. Accessories

5.1. Head Rest

A height and depth adjustable headrest can be mounted on your Eltego.

5.2. Crutch Holder

Put the crutch in the crutch holder and then fix it to the backrest. Never step on the crutch holder and never use it as a tipping aid.

6. Technical Data

Dimensions and Weights

Max. Load	120 kg
Seat width:	400 - 540 mm (via armrest adjustment)
Seat height:	485 mm (without cushion)
Depth of Seat:	440 mm
Seat Angle:	3°
Height of Backrest:	540 mm
Backrest angle:	90° - 120° (in 4 x 10° steps)
Overall length:	1050 mm (including footrests)
Total Width:	600 mm
Empty weight:	72,6 kg (including batteries)
Weight 2 batteries:	21.2 kg

Tires

Drive Wheels:	12,5" x 21/4" PU / 320 x 55 mm PU
Castors:	8" PU / 200 x 50 mm PU

Electrical System

Main fuse:	50 A
Batteries:	2 x 12 V / 36 Ah
Charging Unit:	24 V DC 3 A
Motor:	2 x 24 V / 250 W / 4600 RPM
Controller:	PG Drive VSI
Lighting:	rear lights with separate activation
Lighting (optional):	headlights and rear lights according to the German traffic regulations [StVZO]

Driving Specifications

Speed:	6 km/h (forwards) / 4 km/h (reverse)
Braking distance at 6km/h:	1000 mm*
Max. Incline:	max. 8°

Max. Obstacle Height:	max. 50 mm
Range:	25 km (may be reduced by frequent use on slopes, rough ground, climbing kerbs, etc.)
Turning Radius:	1200 mm

Fittings

Frame:	Powder-coated aluminum frame
Brake:	Automatic braking system, Electromagnetic motor brake according to the German traffic regulations [StVZO] Resistance to ignition according to EN1021-1/2.

Operating forces

Joystick:	approx. 1 N
Buttons on the controller:	approx. 1 N
Disconnecting the drives:	> 60 N

Ambient conditions

Temperature:	-10 – 50°C
Air humidity:	20 – 80%

Storage conditions

Temperature:	-20 – 65°C
Air humidity:	20 – 60%

Accessories

Headrest, cane holder.

*Please note that the braking distance may be longer on uneven, steep surfaces.
The dimensions given are subject to manufacturing tolerances of +/- 10mm.

7. Maintenance Tips for Users

7.1. Cleaning and Disinfection

Cleaning

Make sure that the electric wheelchair is switched off before you start cleaning.

- The frame parts and casings can be wiped with a damp cloth. If there is heavy contamination, you can use a mild household cleaner.
- The wheels can be cleaned with a damp brush with plastic bristles (do not use a wire brush!).
- You can clean dust and minor contamination from your seat, backrest and armrests with a damp cloth. Do not use any aggressive cleaners as these could damage the surface of the material. The upholstery can be washed with 40 °C warm water.
- Make sure that electrical or electronic components do not come into contact with water!

Disinfection

Before disinfection clean the foldable electric wheelchair thoroughly in accordance with the instructions:

- Use a household disinfectant for wipe-down disinfection.
- Spray the upholstery/pads and handles as appropriate.
- Leave the cleaner to work in accordance with the manufacturer's specifications, before you wipe it off.

When doing this, make sure that the electrics and electronics do not come into contact with liquid.

Warning!
Do not use a high-pressure or steam cleaner.

Please observe the enclosed overview on care, hygiene and maintenance work

7.2. Check Before Use

Warning!
Do not use the wheelchair if it is behaving abnormally or erratically.

7.2.1. Brake Inspection

Make sure your brakes are operating properly before each use. Carefully drive off and then release the joystick, the wheelchair must come to a standstill as usual.

When stopped, the mechanical brakes must engage with an audible click. Should you notice any unusual braking behavior, please notify your qualified mobility dealer right away and ensure the wheelchair is not being used.

The motor brake is functioning correctly if you can't push the wheelchair when it is turned off.

Should you nevertheless be able to do so, then the electro-magnetic brakes are defective. In this case you may no longer use the wheelchair. Contact your qualified mobility dealer immediately.

7.2.2. Checking the Lights and Horn

Make sure your lights and horn are functioning properly before each use.

If a unit under inspection does not react when the corresponding button on the control panel is pushed, please notify your qualified mobility dealer immediately and ensure the wheelchair is not being used.

7.2.3. Inspection of Tires

Before embarking on any journey, check that the tires are in good condition. Carry out a regular check of the tread depth of the tires once a month. This should be at least 1 mm.

Any faulty or worn out tires must be changed through the specialist retailer.

7.3. Batteries

Your wheelchair is equipped with maintenance-free gel batteries whose fluid level doesn't need to be checked and which are electronically protected against total discharge.

Should your wheelchair not be used for a longer time, it should be connected to the charging unit every 2 months; otherwise a total discharge could destroy the batteries.

Before driving off, just look at the battery status display on your control panel to check the charging status of your batteries. If only the red LEDs light up in the battery status display, then you shouldn't drive for longer distances but connect the wheelchair to the charging unit immediately.

We recommend connecting the charging unit after every use.

Only use an approved charger to charge the batteries (see technical data). Please also follow the operating instructions for the charger!

7.3.1. Charging Procedure

Only use approved chargers to charge the batteries (type in accordance with the technical data or as delivered by the Bischoff & Bischoff dealer Fig. 27).

Please also follow the operating instructions for the charger!

Do not charge frozen batteries.

The batteries must be charged when the temperature is between 0° and 40°C.

When using the scooter for the very first time, we recommend driving until just 3 bars are visible in the display (see Chapter 6.1) and recharge the batteries at that point.

In the event that this charge level is displayed, you should recharge the batteries as soon as possible in order to avoid deep discharge.

Recharging requires approx. 12 - 16 hours.


When using the product after that, you can recharge the batteries at any battery charge level.

You can interrupt the charging procedure even if the batteries are not fully charged, if you need to continue driving, for example.

But where possible, always charge to 100% so that you have the greatest range.

Charging batteries:

- First switch the scooter off.
- Before you start using the charger, check that the power supply is 230 V.
- Connect the charger with the charging socket of the scooter (Fig. 28) and then to the power socket (230 V).
- The charging socket is located on the side below the control unit (Fig. 29).

- Once the charging procedure is complete, unplug the mains plug and disconnect the scooter from the charger.
- When not in use, the battery should be charged once a month in order to avoid deep discharge.
- Deep discharge shortens the service life of the battery significantly, or the battery may also be damaged by the deep discharge.
- Batteries should be stored upright, in dry, frost-free conditions.

As soon as the range with fully-charged batteries starts to decrease significantly, you should replace the batteries.

7.3.2. Exchanging Batteries

Be extremely careful when working with the batteries. Avoid working with metal tools, do not touch the battery terminals and observe the instructions of the manufacturer.

Warning! This work should only be carried out by a qualified mobility dealer. There is a risk to life!

Only use batteries that are approved for this wheelchair (see Technical Data).

Make sure the wheelchair is turned off before starting to exchange the batteries.

- Loosen the fitting holding the batteries in place and pull out the battery sled.
- Loosen the cover of the battery compartment.
- Loosen the battery cable, first the red one and then the black one. You won't need a tool to do this.
- Carefully lift the batteries out on their straps.
- Lift the new batteries into their designated place and reconnect the terminals. Observe the correct color coding (red and black).

- Secure the battery cover with the Velcro strap and slide the battery compartment back into its original position. It will catch automatically.

When the service life of the batteries has expired, i.e. when they can no longer be charged sufficiently, the batteries may not be disposed of with the ordinary household garbage. Old batteries are taken back by the dealer and disposed of properly when buying a new battery from him!

8. Maintenance Instructions

Typical wear parts are brakes and tyres. Contact your dealer for spare parts.

We recommend that you take your wheelchair to your dealer for inspection at least once a year.

If there are faults or defects, take the wheelchair to your dealer immediately for repair.

The information and documentation required for repair and maintenance shall be supplied to the mobility dealer on request.

If you have questions, you can contact the Bischoff & Bischoff customer service. You can find the contact details for Customer Service on the back cover of this instruction manual.

Check your wheelchair regularly in accordance with the maintenance plan.

9. Maintenance Tips for Retailers

Maintenance must be carried out in accordance with the Bischoff & Bischoff maintenance plan. A maintenance schedule will be provided upon request. A training programme for the specialised trade is also available.

Should you have any questions, please contact the Bischoff & Bischoff customer service.

The wheelchair must be inspected according to the maintenance schedule and must be in proper technical condition.

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11. Disposal

If your wheelchair is no longer needed and you want to dispose it, please contact your authorized mobility dealer.

If you want to dispose of it yourself, please ask your local disposal company for the applicable disposal regulations in your town or city.

12. Use by Others

The Eltego is appropriate for use by others. For use by others, the wheelchair must first be thoroughly cleaned and disinfected according to the accompanying care and hygiene instructions.


The product shall get checked by an authorised mobility dealer for general condition, wear and shall get repaired if necessary.

13. Transport

You can determine whether or not your wheelchair is approved for use as a seat in a vehicle for transporting disabled people by referring to the symbols on the CE sticker on the product:

Wheelchair suitable for use as a seat in a vehicle for transporting disabled people (fig. 30)

Wheelchair not suitable for use as a seat in a vehicle for transporting disabled people (fig. 31)

14. Storage and Shipping

If the wheelchair is to be stored or shipped, all loose parts which are not fixed, except for the drive wheels, must be removed and packed individually in an appropriate cardboard box. The power circuit must be open when stored for long periods of time. The individually packed parts can then be packaged together in a larger cardboard box. We recommend keeping the original packaging and storing it with this purpose in mind, so that it is available if required. Your wheelchair then has optimum protection against environmental influences during storage or transport.

For storage of your commode chair, please observe the information in the technical data. To avoid mould growth and damage to the the

upholstery parts, the storage location should be as dry as possible and not be exposed to direct sunlight.

15. Service life

The expected service life of the Ejoy RD is up to 5 years, if it is used correctly, and safely, following the care and maintenance instructions. Beyond this period the wheelchair can be used if it is in a condition allowing its safe usage.

16. Warranty

1. For products supplied by Bischoff & Bischoff, the period of warranty is 24 months from purchase date. The products are delivered free from production and material defects. In case of evidence of material or manufacturing defect, defective parts will be replaced free of charge.

2. This warranty excludes claims, if a repair or replacement of a product or a part is required for the following reasons:

- a) normal wear and tear, which include but is not limited to the following parts where fitted; Batteries, motor brushes, handles, armrests, upholstery, tyres, brakes, caps, etc.
- b) overloading the product such as, for example, by exceeding the maximum user weight or load.
- c) the product or part has not been maintained or serviced in accordance with the Instruction manual, the care and hygiene instructions or the manufacturer's recommendations as shown in the maintenance instructions.
- d) Accessories have been used which are not original accessories.
- e) The product or part has been damaged through negligence, accident or incorrect use.
- f) changes/modifications have been made to the product or parts, which deviate from the manufacturer's specifications.
- g) Repairs have been carried out, before our customer service has been informed of the circumstances.

3. To enforce the warranty, please contact your dealer with the exact details of the nature of the difficulty. The product must be repaired by

a dealer authorised by Bischoff & Bischoff.

4. For parts, which have been repaired or exchanged within the scope of the warranty, the warranty will be extended to the remaining warranty period for the product in accordance with point 1).

5. For original spare parts which have been fitted at the customer's expense, these will have a 12 months guarantee, following the fitting.

6. If the repair fails after a reasonable time, the customer can, at his discretion, request a reduction, compensation or the termination of the purchase contract.

7. The guarantee is subject to the law of the country in which the product was purchased from Bischoff & Bischoff & Bischoff.

Please also observe the Bischoff & Bischoff care, warranty, hygiene and maintenance instructions. Your dealer will be happy to provide these to you.

For replacement devices and repairs, the warranty period is one year.

Warning! Non-observance of the instruction manual, improper maintenance and in particular any technical changes and supplements (added fittings) without the approval of Bischoff & Bischoff will generally lead to forfeiture of warranty and product liability claims.

Warranty-Certificate

(Please note our general business terms and conditions)

Model: Eltego

Serial Number:*

* (to be filled in by authorized dealer)

Dealer:

Date and Stamp

*The serial number is located on the side frame

Note!

In spite of careful research and editing, incorrect information may have crept into our operating manuals. Rights to technical changes reserved. Illustrations may include accessories. The current version of each instruction manual is ready for download at the Bischoff & Bischoff website:
www.bischoff-bischoff.com

Notes

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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